

Kodak

CleAR+ Lens




Help Your Patients See Life Clearly.

Engineered for Performance and Protection

Both sides of KODAK CleAR+ Lens include:

- Protective layers that provide durability and scratch resistance.
- NEW!** Anti-static properties that resist clinging dust.
- Several microscopic layers of minerals that effectively reduce reflections over a broad band of wavelengths.
- Hydrophobic top coat seals the lens to help resist dirt and smudges.

With ordinary lenses, approximately 10% of light is reflected away and doesn't reach the eyes. More reflections and less light mean poorer vision. With KODAK CleAR+ Lens, approximately 99% of available light reaches the eyes providing clearer and crisper vision.



Regular Lens



KODAK CleAR+ Lens

Kodak CleAR+ Lens

Easy to Dispense

KODAK CleAR+ Lens allows eyecare professionals to offer patients a recognizable brand known for trust, respect and loyalty.

The brand is KODAK.
The quality you and
your patients trust.

KODAK CleAR+ Lens carries one of the most respected brand names in the world. KODAK CleAR+ Lens was designed and tested by Signet Armorlite, a company with over 66 years of leadership and innovation in ophthalmic lens technology.



Easy-to-Clean

KODAK CleAR+ Lens resists oily smudges, dust and dirt, making the lens easy to clean.



Scratch Resistant

KODAK CleAR+ Lens is protected from irritating scratches so the lens lasts longer and looks better over the life of the prescription.



Eliminates Irritating Glare

KODAK CleAR+ Lens reduces eyestrain and fatigue from artificial lighting and computer screens.




Reduces Nighttime Reflections

KODAK CleAR+ Lens improves ability to see objects in low light for safer, more comfortable night driving.

See the *Colors* of Life

Kodak LENS



Signet Armorlite
5803 Newton Drive, Suite A
Carlsbad, CA 92008
800-759-0075
www.signetarmorlite.com



Photographic demonstrations simulated.
The Kodak trademark and trade dress are used under license from Kodak by Signet Armorlite, Inc. CleAR is a trademark of Signet Armorlite, Inc.
© 2015 Signet Armorlite, Inc.