



Navigator

Navigator® and
Navigator Short® Progressives

Great performance
at an excellent price.

Soft Design Progressive

Easy Adaptation

Two Corridor Lengths

Navigator Progressives are created with an asymmetrical design to deliver a wide distance viewing area and broad corridor ideal for first time progressive users. At a great value, it is also an excellent option for insurance jobs and tight budgets. Navigator lenses are available in regular and short corridor designs to accommodate most frame styles.

PhotoViews™ is a photochromic lens option that provides fast fade time from dark to clear and 100% UV protection on the front surface of the lens.

Lens Material	Diameter	Base Curves	Add Powers	Power Range
Navigator (20mm minimum fitting height)				
Signia™ 1.50 (coated and uncoated)	78mm	2.00, 4.50, 5.75, 7.75	.75 to 3.50 in .25D steps	-8.50 to +5.50
PhotoViews 1.50	76mm	2.00, 4.50, 5.75, 7.75	1.00 to 3.00 in .25D steps	-8.00 to +5.00
Signia Polarized 1.50	78mm	3.25, 5.00, 7.00	1.00 to 3.00 in .25D steps	-7.00 to +5.00
Signia Poly	78mm	2.00, 4.50, 5.75, 7.75	1.00 to 3.00 in .25D steps	-10.00 to +6.00
Signia Lite 1.60	78mm	.50, 2.00, 4.50, 5.75, 7.75	1.00 to 3.00 in .25D steps	-11.75 to +6.50
Navigator Short (14mm minimum fitting height)				
Signia™ 1.50 (coated and uncoated)	74mm	1.25, 3.25, 5.50, 7.50	1.00 to 3.00 in .25D steps	-9.00 to +5.00
PhotoViews 1.50	74mm	1.25, 3.25, 5.50, 7.50	1.00 to 3.00 in .25D steps	-8.00 to +5.00
Signia Poly	78mm	2.00, 4.50, 5.75, 7.25	1.00 to 3.00 in .25D steps	-10.00 to +6.00

Note: Signia and PhotoViews materials are compatible with all major AR coatings.

Signet Armorlite • 5803 Newton Drive, Suite A, Carlsbad, CA 92008 • 800-759-0075 • www.signetarmorlite.com